

new year solutions

£15,000 & Under



What is it? Phaser-EPL

What does it do? Treats unwanted hair, even fine and lighter coloured hairs, and rejuvenates skin including unwanted pigmentation.

How does it work? Phaser-EPL is an Intense Pulsed Light System utilising its own Smartlite™ Dual Phase Technology, which produces high levels of light energy. It uses a combination of enclosed waveguide and a special reflecting cavity at the rear of the flash lamp source which increases transmission of light energy into tissue significantly compared to a standard IPL. The Dual Phase Pulse destructs both the bulge and the bulb for ultimate hair removal results.

Trade price: On application (lease payments available)

• Chromogenex Technologies Ltd:
01554 755 444

£20,000 & Under



What is it? Lynton Excelight

What does it do? Delivers advanced Intense Pulsed Light (IPL) treatments including hair removal, skin rejuvenation, anti-ageing, Acne, red vein removal, sun damage, pigmentation, rosacea and more.

How does it work? Performing over 10 different treatment areas, this is Lynton's most advanced multi-treatment IPL system. It features unique, ergonomic ILG™ (Interchangeable Light Guide) technology which gives you flexibility in treatment spot size, with a 'snap in' design to ensure quick light guide replacement, and quick, comfortable treatments.

Trade price: On application

• Lynton: 01477 536 977



What is it? eSpa™

What does it do? Tightens skin on the face and body.

How does it work? This dual-application platform is powered by proprietary eLos™ technology; a combination of gentle pulsed infrared light (wavelengths 700-2000nm) and Bi-polar radio frequency (RF) provide deep dermal heating, up to 3 mm in depth. New collagen is formed via a wound healing response and visible results are revealed immediately.

It features the Motif Mode with fast energy delivery; the low fluence and high repetition rate enable an efficient pain free and safe treatment for all skin types.

Trade price: £17,995 + VAT

• Syneron Candela: 0845 521 0698

What is it? Century 808 Diode Laser

What does it do? Standard and fast Laser hair removal on all clients including those with darker skins and thick darker hair, as well as blonde hair on fairer skins.

How does it work? The Century uses an optimum wavelength to destroy the hair follicle – a long pulse 808nm frequency – with a high repetition rate allowing for a fast treatment time with the added benefit of an adjustable cooling head. Clients are said to experience permanent hair reduction in 100 days (five sessions, once every 20 days).

Trade price: £16,500 + VAT

• Premier Lasers: 0800 043 7272



* Prices are for guidance only and are subject to change.
Finance options subject to status.
Prices may or may not include VAT.