

## PlumeSafe® Turbo - Surgical Smoke Evacuator

PlumeSafe® Turbo is a powerful portable surgical smoke evacuator designed with sufficient flow for the aesthetic market. It features turbo suction power, an interactive LCD display and filter tracking technology. It can easily be placed on a counter, shelf, floor or other location convenient to the operative area.

100/120 Volt 50/60 Hz .....PSWTURBO AS  
220/240 Volt 50/60 Hz .....PSWTURBO\_AS/220  
(includes one footswitch; filter sold separately)

PlumeSafe® Turbo Replacement Filters  
VS06, VSXLL, VSXLLT, VSEY, VSOH, and VSVET Series  
Filter Life ..... 6-24 hours (Depending on Filter model and flow rate)



Pneumatic Footswitch  
*(included)*